



INDEPENDENCE SERIES 3-WHEEL SCOOTER



- ✓ Weighs 66 Pounds w/o battery & can be separated into two halves (Front Chassis 35 Pounds, Rear Chassis 31 Pounds)
- ✓ Folds to the size of a carry on Suitcase (29" x 19" x 15")
- ✓ Supports Weight of up to 265 Pounds
- ✓ Travels up to 24 Miles
- ✓ Speed up to 3.75 Miles Per Hour (6 speeds)
- ✓ 250 Watt Brushless DC Motor
- ✓ Intelligent Electromagnetic Brakes
- ✓ 48V/20ah Lithium Ion Battery
- ✓ Includes USB Charger, Mobile Phone Holder, LED Light
- ✓ 2 Year Limited Warranty



250 Watt Brushless DC Motors



6 Speed Adjustments



Steering Column Height Adjustments



USB Charger



Steering Column Angle Adjustments



Horn



Solid Polyurethane/Rubber 8" And 9" Tires



Battery Power Indicator Light



Forward/Reverse Controls



Light



Removable Lithium Ion Battery 48V/20ah



Folds To The Size Of A Suitcase



Separates Into 2 Halves For Easy Lifting



Wide Footpad



3 Seat Height Adjustments



Weight Capacity 265 lbs



Scooter Weight (w/o battery) 66 lbs



Speed 3.75 mph



Travel Distance 25 Miles



Power 250 Watts

INDEPENDENCE SERIES

Miracle Mobility

Miracle Mobility™



Weight Capacity
265 lbs



Scooter Weight
(w/o battery)



Speed
3.75 mph



Travel Distance
25 Miles



Power
250 Watts

Miracle Mobility Independence 3-Wheel Scooter 2019 Model (discontinued) (White)

Color Options	White
Model #	MM-IS3WS-W3 (X1S)
UPC Code	7-02444-00734-7 (white)
Motor Type	Magnetic Brushless DC Motor
Motor Power Rating	250 Watts
Weight Capacity	265 lbs
Chair Weight w/o Battery	66 lbs
Chair Weight w/Battery	72 lbs
Weight Separated	35 lbs (front), 31 lbs (back)
Maximum Speed	3.75 mph
Speed Settings	6 Speeds
Driving Distance	up to 25 miles
Incline Grade Capacity	10 degrees
Battery Type	Lithium Ion Battery
Battery Capacity	48V/20AH
Batteries Per Unit	1 Batteries
Charge Time	4 hours
Joystick Controller	DC24V/2A
Charger	24V/2A
Brakes	Intelligent Electromagnetic
Seat Size	14" x 14.5"
Seat Height	12.5" H
Tire Material	Solid Polyurethane & Rubber
Front Tire Size	8"
Rear Tire Size	9"
Frame Construction	Aluminum
Gross Shipping Weight	75 lbs
Net Weight	66 lbs
Item Dimensions	33" H x 25" W x 48" L
Item Dimensions folded	29" H x 19" W x 15" L
Carton Dimensions	21" x 17" x 32"
Carton Cube	6.6
Warranty	2 yrs (frame), 1 yr (motor), 1 yr (parts)
FDA Approved	no
Insurance Code	NA