


Miracle Mobility Freedom Series Gold (D05) 2019 Model

Image	
Color Options	Silver
Model #	MM-FSG-W1 (silver)
UPC Code	7-02444-00758-3 (silver)
Motor Type	Magnetic Brushless DC Motor
Motor Power Rating	180 Watts x 2 (360 Watts Total)
Weight Capacity	265 lbs
Chair Weight w/o Battery	41 lbs
Chair Weight w/Battery	46 lbs
Maximum Speed	3.75 mph
Speed Settings	5 speeds
Driving Distance	up to 12 miles
Incline Grade Capacity	12 degrees
Battery Type	Lithium Ion Battery
Battery Capacity	24V/10AH/1pcs
Batteries Per Unit	1 Battery
Charge Time	6-8 hours
Joystick Controller	DC24V/2A
Charger	24V/2A
Brakes	Intelligent Electromagnetic
Seat Size	15.75" D x 15.75" W
Seat Height	19" H
Tire Material	Solid Polyurethane & Rubber
Front Tire Size	7"
Rear Tire Size	8"
Frame Construction	Aluminum
Gross Shipping Weight	60 lbs
Net Weight	48 lbs
Item Dimensions	40" L x 36" H x 23" W
Item Dimensions folded	13" L x 31" H x 24" W
Carton Dimensions	26" x 17" x 31"
Carton Cube	7.9
Warranty	5 yrs (frame), 2 yrs (motor), 1 yr (parts)
FDA Approved	Yes
Insurance Code	K0824, K0813