


Miracle Mobility Freedom Series Gold (D20) New 2020 Model

Image	
Color Options	Black & White Pink & Silver Black & Silver
Model #	MM-FSG-W9D20 (Black & White) MM-FSG-W8D20 (Pink & Silver) MM-FSG-W10D20 (Black & Silver)
UPC Code	8-50012-19833-8 (Black & White) 8-50012-19835-2 (Pink & Silver) 8-50012-19834-5 (Black & Silver)
Motor Type	Magnetic Brushless DC Motor
Motor Power Rating	180 Watts x 2 (360 Watts Total)
Weight Capacity	265 lbs
Chair Weight w/o Battery	38 lbs
Chair Weight w/Battery	44 lbs
Maximum Speed	3.75 mph
Speed Settings	5 speeds
Driving Distance	up to 12 miles
Incline Grade Capacity	12 degrees
Battery Type	Lithium Ion Battery
Battery Capacity	24V/10AH/1pcs
Batteries Per Unit	1 Battery
Charge Time	6-8 hours
Joystick Controller	DC24V/2A
Charger	24V/2A
Brakes	Intelligent Electromagnetic
Seat Size	16" D x 16" W
Seat Height	19" H
Tire Material	Solid Polyurethane & Rubber
Front Tire Size	7"
Rear Tire Size	8"
Frame Construction	Aluminum
Gross Shipping Weight	58 lbs
Net Weight	46 lbs
Item Dimensions	36" L x 36" H x 23" W
Item Dimensions folded	12" L x 28" H x 25" W
Carton Dimensions	29" x 33" x 16.5"
Carton Cube	9.1
Warranty	5 yrs (frame), 2 yrs (motor), 1 yr (parts)
FDA Approved	YES
Insurance Code	K0824, K0813