

Miracle Mobility Freedom Series Gold XL D10 New 2020 Model



Color Options	Silver - W1D10 Black & White - W9D10 Black & Blue - W5D10 Silver & Pink - W8D10
Model #	MM-FSGXL-W1D10 (Silver) MM-FSGXL-W9D10 (Black & White) MM-FSGXL-W5D10 (Black & Blue) MM-FSGXL-W8D10 (Pink & Silver)
UPC Code	8-50012-19802-4 (Silver) 8-50012-19831-4 (Black & White) 8-50012-19832-1 (Black & Blue) 8-50012-19830-7 (Gunmetal & Pink)
Motor Type	Magnetic Brushless DC Motor
Motor Power Rating	250 Watts x 2 (500 Watts Total)
Weight Capacity	396 LBS
Chair Weight w/o Battery	55 lbs
Chair Weight w/Battery	60 lbs
Maximum Speed	3.75 mph
Speed Settings	5 Speeds
Driving Distance	14 miles (2 batteries) 25 miles (4 batteries)
Incline Grade Capacity	10 degrees
Battery Type	Lithium Ion Battery
Battery Capacity	24V/6AH/2pcs
Batteries Per Unit	2 Batteries
Battery Housing Case	Up to 4 Batteries
Charge Time	4-6 hours
Joystick Controller	DC24V/2A
Charger	24V/2A
Brakes	Intelligent Electromagnetic
Seat Size	17" D x 17" W
Seat Height	20" H
Tire Material	Solid Polyurethane & Rubber
Front Tire Size	8"
Rear Tire Size	12"
Frame Construction	Aluminum Alloy
Gross Shipping Weight	72 lbs
Net Weight	62 lbs
Item Dimensions	37" L x 38" H x 23" W
Item Dimensions folded	18" L x 28" H x 24" W
Carton Dimensions	29" x 34" x 18"
Carton Cube	10.1
Warranty	5 yrs (frame), 2 yrs (motor), 1 yr (parts)
FDA Approved	Yes
Insurance Code	K0824, K0813