


## Miracle Mobility Independence 3-Wheel Scooter 2019 Model (White)

|                          |  |
|--------------------------|--|
| <b>Image</b>             |  |
| Color Options            | White  |
| Model #                  | MM-IS3WS-W3 (X1S)  |
| UPC Code                 | 7-02444-00734-7 (white)  |
| Motor Type               | Magnetic Brushless DC Motor  |
| Motor Power Rating       | 250 Watts  |
| Weight Capacity          | 265 lbs  |
| Chair Weight w/o Battery | 62 lbs   |
| Chair Weight w/Battery   | 67 lbs   |
| Weight Separated         | 35 lbs (front), 32 lbs (back)  |
| Maximum Speed            | 3.75 mph   |
| Speed Settings           | 6 Speeds   |
| Driving Distance         | up to 25 miles   |
| Incline Grade Capacity   | 10 degrees   |
| Battery Type             | Lithium Ion Battery  |
| Battery Capacity         | 48V/8.8AH  |
| Batteries Per Unit       | 1 Batteries  |
| Charge Time              | 4 hours  |
| Joystick Controller      | DC24V/2A   |
| Charger                  | 24V/2A   |
| Brakes                   | Intelligent Electromagnetic  |
| Seat Size                | 14" x 14.5"  |
| Seat Height              | 12.5" H  |
| Tire Material            | Solid Polyurethane & Rubber  |
| Front Tire Size          | 8"   |
| Rear Tire Size           | 9"   |
| Frame Construction       | Aluminum   |
| Gross Shipping Weight    | 75 lbs   |
| Net Weight               | 67 lbs   |
| Item Dimensions          | 33" H x 25" W x 48" L  |
| Item Dimensions folded   | 29" H x 20" W x 12" L  |
| Carton Dimensions        | 21" x 17" x 32"  |
| Carton Cube              | 6.6  |
| Warranty                 | 2 yrs (frame), 1 yr (motor), 1 yr (parts)  |
| FDA Approved             | no   |
| Insurance Code           | NA   |